

SERPING1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:1000 - 1:2000 IHC 1:50 - 1:200

Calculated MW:55kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human SERPING1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C1IN; C1NH; HAE1; HAE2; C1INH;

Background:

This gene encodes a highly glycosylated plasma protein involved in the regulation of the complement cascade. Its protein inhibits activated C1r and C1s of the first complement component and thus regulates complement activation. Deficiency of this protein is associated with hereditary angioneurotic oedema (HANE). Alternative splicing results in multiple transcript variants encoding the same isoform. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).

Catalog #:A1717

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:710

Isotype:IgG

Swiss Prot:P05155

Purity:Affinity purification

For research use only.