

## IGFBP5

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:100

**Calculated MW:**31kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human IGFBP5

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

IGFBP5; IBP5;

**Background:**

IGFBP5, also named as IBP5, prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. It alter the interaction of IGFs with their cell surface receptors. This antibody is specific to IGFBP5.

*To place an order, please [Click HERE](#).*

**Catalog #:**A1720

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**3488

**Isotype:**IgG

**Swiss Prot:**P24593

**Purity:**Affinity purification

For research use only.