

PGF

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:24kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PGF

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

a

Synonym:

PGFL; PLGF; PIGF-2; D12S1900; SHGC-10760;

Catalog #:A1727

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5228

Isotype:IgG

Swiss Prot:P49763

Purity:Affinity purification

For research use only.

Background:

This gene encodes a growth factor found in placenta which is homologous to vascular endothelial growth factor. Alternatively spliced transcripts encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).