

SCGB1A1

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:10kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human SCGB1A1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UGB; UP1; CC10; CC16; CCSP; UP-1; CCPBP;

Background:

This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma. [provided by RefSeq, May 2010]

To place an order, please [Click HERE](#).

Catalog #:A1730

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7356

Isotype:IgG

Swiss Prot:P11684

Purity:Affinity purification

For research use only.