

PTH1R

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:66kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PTH1R

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PFE; PTHR; PTHR1;

Background:

The protein encoded by this gene is a member of the G-protein coupled receptor family 2. This protein is a receptor for parathyroid hormone (PTH) and for parathyroid hormone-like hormone (PTHrP). The activity of this receptor is mediated by G proteins which activate adenylyl cyclase and also a phosphatidylinositol-calcium second messenger system. Defects in this receptor are known to be the cause of Jansen's metaphyseal chondrodysplasia (JMC), chondrodysplasia Blomstrand type (BOCD), as well as enchondromatosis. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, May 2010]

To place an order, please [Click HERE](#).

Catalog #:A1744

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5745

Isotype:IgG

Swiss Prot:Q03431

Purity:Affinity purification

For research use only.