

MTTP

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 99kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MTTP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ABL; MTP;

Catalog #: A1746

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4547

Isotype: IgG

Swiss Prot: P55157

Purity: Affinity purification

For research use only.

Background:

MTP encodes the large subunit of the heterodimeric microsomal triglyceride transfer protein. Protein disulfide isomerase (PDI) completes the heterodimeric microsomal triglyceride transfer protein, which has been shown to play a central role in lipoprotein assembly. Mutations in MTP can cause abetalipoproteinemia. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).