

## FKBP1A

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**12kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human FKBP1A

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

FKBP1; PKC12; PKCI2; FKBP12; PPIASE; FKBP-12; FKBP-1A;

**Background:**

The protein encoded by this gene is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. The protein is a cis-trans prolyl isomerase that binds the immunosuppressants FK506 and rapamycin. It interacts with several intracellular signal transduction proteins including type I TGF-beta receptor. It also interacts with multiple intracellular calcium release channels, and coordinates multi-protein complex formation of the tetrameric skeletal muscle ryanodine receptor. In mouse, deletion of this homologous gene causes congenital heart disorder known as noncompaction of left ventricular myocardium. Multiple alternatively spliced variants, encoding the same protein, have been identified. The human genome contains five pseudogenes related to this gene, at least one of which is transcribed. [provided by RefSeq, Sep 2008]

*To place an order, please [Click HERE](#).*

**Catalog #:**A1763

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**2280

**Isotype:**IgG

**Swiss Prot:**P62942

**Purity:**Affinity purification

For research use only.