

HGS

Reactivity:Human Mouse

Tested applications:WB IF IP

Recommended Dilution:WB 1:500 - 1:2000 IF 1:50 - 1:200 IP 1:20 - 1:50

Calculated MW:86kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human HGS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HRS;

Background:

The protein encoded by this gene regulates endosomal sorting and plays a critical role in the recycling and degradation of membrane receptors. The encoded protein sorts monoubiquitinated membrane proteins into the multivesicular body, targeting these proteins for lysosome-dependent degradation.

To place an order, please [Click HERE](#).

Catalog #:A1790

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9146

Isotype:IgG

Swiss Prot:O14964

Purity:Affinity purification

For research use only.