

DCK

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:31kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human DCK

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MGC117410; MGC138632;

Background:

Deoxycytidine kinase (dCK) is required for the phosphorylation of several deoxyribonucleoside analogues that are widely employed as chemotherapeutic agents(PMID:12628850). It is also a key enzyme in the phosphorylation of a variety of antineoplastic and antiviral nucleoside analogs including cytosine arabinoside (araC) and dideoxycytidine (ddCyd); deficiency of deoxycytidine kinase activity mediates resistance to these drugs(PMID:1996353). This protein can exist as a homodimer(PMID:16421443). This antibody is specific to DCK.

To place an order, please [Click HERE](#).

Catalog #:A1794

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1633

Isotype:IgG

Swiss Prot:P27707

Purity:Affinity purification

For research use only.