

## CXADR

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**40kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human CXADR

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CAR; HCAR; CAR4/6;

**Background:**

The protein encoded by this gene is a type I membrane receptor for group B coxsackieviruses and subgroup C adenoviruses. Several transcript variants encoding different isoforms have been found for this gene. Pseudogenes of this gene are found on chromosomes 15, 18, and 21.

*To place an order, please [Click HERE](#).*

**Catalog #:**A1822

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**1525

**Isotype:**IgG

**Swiss Prot:**P78310

**Purity:**Affinity purification

For research use only.