

FCER2

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:36kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human FCER2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD23; FCE2; CD23A; IGEBF; CLEC4J;

Background:

The protein encoded by this gene is a B-cell specific antigen, and a low-affinity receptor for IgE. It has essential roles in B cell growth and differentiation, and the regulation of IgE production. This protein also exists as a soluble secreted form, then functioning as a potent mitogenic growth factor. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.[provided by RefSeq, Jul 2011]

To place an order, please [Click HERE](#).

Catalog #:A1839

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2208

Isotype:IgG

Swiss Prot:P06734

Purity:Affinity purification

For research use only.