

POLR2J

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:13kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human POLR2J

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

ei

Synonym:

RPB11; RPB11A; RPB11m; hRPB14; POLR2J1;

Catalog #:A1843

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5439

Isotype:IgG

Swiss Prot:P52435

Purity:Affinity purification

For research use only.

Background:

This gene encodes a subunit of RNA polymerase II, the polymerase responsible for synthesizing messenger RNA in eukaryotes. The product of this gene exists as a heterodimer with another polymerase subunit; together they form a core subassembly unit of the polymerase. Two similar genes are located nearby on chromosome 7q22.1 and a pseudogene is found on chromosome 7p13.

To place an order, please [Click HERE](#).