

RNASE3

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:200 IHC 1:50 - 1:100

Calculated MW:18kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RNASE3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ECP; RNS3;

Background:

The protein encoded by this gene belongs to the pancreatic ribonuclease family, a subset of the ribonuclease A superfamily. The protein exhibits antimicrobial activity against pathogenic bacteria

To place an order, please [Click HERE](#).

Catalog #:A1854

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6037

Isotype:IgG

Swiss Prot:P12724

Purity:Affinity purification

For research use only.