

POLR2I

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:15kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human POLR2I

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

ei

Synonym:

RPB9; hRPB14.5;

Catalog #:A1865

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5438

Isotype:IgG

Swiss Prot:P36954

Purity:Affinity purification

For research use only.

Background:

This gene encodes a subunit of RNA polymerase II, the polymerase responsible for synthesizing messenger RNA in eukaryotes. This subunit, in combination with two other polymerase subunits, forms the DNA binding domain of the polymerase, a groove in which the DNA template is transcribed into RNA. The product of this gene has two zinc finger motifs with conserved cysteines and the subunit does possess zinc binding activity. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).