

GABRB2

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:59kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human GABRB2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

g

Catalog #:A1876

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2561

Isotype:IgG

Swiss Prot:P47870

Purity:Affinity purification

For research use only.

Background:

The gamma-aminobutyric acid (GABA) A receptor is a multisubunit chloride channel that mediates the fastest inhibitory synaptic transmission in the central nervous system. This gene encodes GABA A receptor, beta 2 subunit. It is mapped to chromosome 5q34 in a cluster comprised of genes encoding alpha 1 and gamma 2 subunits of the GABA A receptor. Alternative splicing of this gene generates 2 transcript variants, differing by a 114 bp insertion.

To place an order, please [Click HERE](#).