

RXRG

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 51kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RXRG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RXRC; NR2B3;

Catalog #: A1877

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6258

Isotype: IgG

Swiss Prot: P48443

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the retinoid X receptor (RXR) family of nuclear receptors which are involved in mediating the antiproliferative effects of retinoic acid (RA). This receptor forms dimers with the retinoic acid, thyroid hormone, and vitamin D receptors, increasing both DNA binding and transcriptional function on their respective response elements. This gene is expressed at significantly lower levels in non-small cell lung cancer cells. Alternatively spliced transcript variants have been described.

To place an order, please [Click HERE](#).