

RXRG

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:51kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RXRG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RXRC; NR2B3;

Catalog #:A1877

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6258

Isotype:IgG

Swiss Prot:P48443

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the retinoid X receptor (RXR) family of nuclear receptors which are involved in mediating the antiproliferative effects of retinoic acid (RA). This receptor forms dimers with the retinoic acid, thyroid hormone, and vitamin D receptors, increasing both DNA binding and transcriptional function on their respective response elements. This gene is expressed at significantly lower levels in non-small cell lung cancer cells. Alternatively spliced transcript variants have been described.

To place an order, please [Click HERE](#).