

ZFP36

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:34kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human ZFP36

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

TTP; G0S24; GOS24; TIS11; NUP475; zfp-36; RNF162A;

To place an order, please [Click HERE](#).

Catalog #:A1878

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7538

Isotype:IgG

Swiss Prot:P26651

Purity:Affinity purification

For research use only.