

EPHX2

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW:62kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human EPHX2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CEH; SHE; Epoxide Hydratase;

Background:

This gene encodes a member of the epoxide hydrolase family. The protein, found in both the cytosol and peroxisomes, binds to specific epoxides and converts them to the corresponding dihydrodiols. Mutations in this gene have been associated with familial hypercholesterolemia. Alternatively spliced transcript variants have been described.

To place an order, please [Click HERE](#).

Catalog #:A1885

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2053

Isotype:IgG

Swiss Prot:P34913

Purity:Affinity purification

For research use only.