

## EPHX2

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:1000 IHC 1:50 - 1:100

**Calculated MW:** 62kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human EPX2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CEH; SHE; Epoxide Hydratase;

**Catalog #:** A1885

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 2053

**Isotype:** IgG

**Swiss Prot:** P34913

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a member of the epoxide hydrolase family. The protein, found in both the cytosol and peroxisomes, binds to specific epoxides and converts them to the corresponding dihydrodiols. Mutations in this gene have been associated with familial hypercholesterolemia. Alternatively spliced transcript variants have been described.

*To place an order, please [Click HERE](#).*