

LIPG

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 57kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LIPG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

EL; EDL; PRO719;

Catalog #: A1891

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9388

Isotype: IgG

Swiss Prot: Q9Y5X9

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene has substantial phospholipase activity and may be involved in lipoprotein metabolism and vascular biology. This protein is designated a member of the TG lipase family by its sequence and characteristic lid region which provides substrate specificity for enzymes of the TG lipase family.

To place an order, please [Click HERE](#).