

SRD5A1

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:29kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SRD5A1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

S5AR 1;

Catalog #:A1898

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6715

Isotype:IgG

Swiss Prot:P18405

Purity:Affinity purification

For research use only.

Background:

Steroid 5-alpha-reductase (EC 1.3.99.5) catalyzes the conversion of testosterone into the more potent androgen, dihydrotestosterone (DHT). Also see SRD5A2 (MIM 607306).

To place an order, please [Click HERE](#).