

PADI4

Reactivity:Human Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:50

Calculated MW:74kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PADI4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

PAD; PAD4; PDI4; PDI5; PADI5;

Catalog #:A1906

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23569

Isotype:IgG

Swiss Prot:Q9UM07

Purity:Affinity purification

For research use only.

Background:

This gene is a member of a gene family which encodes enzymes responsible for the conversion of arginine residues to citrulline residues. This gene may play a role in granulocyte and macrophage development leading to inflammation and immune response.

To place an order, please [Click HERE](#).