

BST2

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:19kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human BST2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD317; TETHERIN;

Catalog #:A1914

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:684

Isotype:IgG

Swiss Prot:Q10589

Purity:Affinity purification

For research use only.

Background:

Bone marrow stromal cells are involved in the growth and development of B-cells. The specific function of the protein encoded by the bone marrow stromal cell antigen 2 is undetermined; however, this protein may play a role in pre-B-cell growth and in rheumatoid arthritis.

To place an order, please [Click HERE](#).