

MNAT1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:36kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human MNAT1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

el

Synonym:

MAT1; TFB3; CAP35; RNF66;

Catalog #:A1918

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4331

Isotype:IgG

Swiss Prot:P51948

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene, along with cyclin H and CDK7, forms the CDK-activating kinase (CAK) enzymatic complex. This complex activates several cyclin-associated kinases and can also associate with TFIIF to activate transcription by RNA polymerase II. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011]

To place an order, please [Click HERE](#).