

## EBAG9

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 24kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human EBAG9

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

b

**Synonym:**

EB9; PDAF;

**Catalog #:** A1935

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 9166

**Isotype:** IgG

**Swiss Prot:** O00559

**Purity:** Affinity purification

For research use only.

**Background:**

This gene was identified as an estrogen-responsive gene. Regulation of transcription by estrogen is mediated by estrogen receptor, which binds to the estrogen-responsive element found in the 5'-flanking region of this gene. The encoded protein is a tumor-associated antigen that is expressed at high frequency in a variety of cancers. Alternate splicing results in multiple transcript variants. A pseudogene of this gene has been defined on chromosome 10.

**To place an order, please [Click HERE](#).**