

## PRDM1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**92kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human PRDM1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

BLIMP1; PRDI-BF1;

**Catalog #:**A1960

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**639

**Isotype:**IgG

**Swiss Prot:**O75626

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a protein that acts as a repressor of beta-interferon gene expression. The protein binds specifically to the PRDI (positive regulatory domain I element) of the beta-IFN gene promoter. Transcription of this gene increases upon virus induction. Two alternatively spliced transcript variants that encode different isoforms have been reported.

*To place an order, please [Click HERE](#).*