

## MMP8

---

**Reactivity:**Human Mouse

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:**53kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human MMP8

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A1963

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**4317

**Isotype:**IgG

**Swiss Prot:**P22894

**Purity:**Affinity purification

For research use only.