

NDP

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:15kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human NDP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ND; EVR2; FEVR

Catalog #:A1964

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4693

Isotype:IgG

Swiss Prot:Q00604

Purity:Affinity purification

For research use only.

Background:

This gene encodes a secreted protein with a cystein-knot motif that activates the Wnt/beta-catenin pathway. The protein forms disulfide-linked oligomers in the extracellular matrix. Mutations in this gene result in Norrie disease and X-linked exudative vitreoretinopathy.

To place an order, please [Click HERE](#).