

## EXT2

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 82kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human EXT2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

SOTV;

**Catalog #:** A1973

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 2132

**Isotype:** IgG

**Swiss Prot:** Q93063

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes one of two glycosyltransferases involved in the chain elongation step of heparan sulfate biosynthesis. Mutations in this gene cause the type II form of multiple exostoses. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

*To place an order, please [Click HERE](#).*