

EXT2

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:82kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human EXT2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SOTV;

Catalog #:A1973

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2132

Isotype:IgG

Swiss Prot:Q93063

Purity:Affinity purification

For research use only.

Background:

This gene encodes one of two glycosyltransferases involved in the chain elongation step of heparan sulfate biosynthesis. Mutations in this gene cause the type II form of multiple exostoses. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

To place an order, please [Click HERE](#).