

HSD17B2

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:42kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human HSD17B2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HSD17; SDR9C2; EDH17B2;

To place an order, please [Click HERE](#).

Catalog #:A1983

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3294

Isotype:IgG

Swiss Prot:P37059

Purity:Affinity purification

For research use only.