

GAL

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:13kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GAL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GALN; GLNN; GMAP; GAL-GMAP;

Background:

Galanin is small neuropeptide that functions as a cellular messenger within the central and peripheral nervous systems, modulating diverse physiologic functions (Mechenthaler, 2008

To place an order, please [Click HERE](#).

Catalog #:A1991

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:51083

Isotype:IgG

Swiss Prot:P22466

Purity:Affinity purification

For research use only.