

CD244

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:41kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CD244

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

2B4; NAIL; Nmrk; NKR2B4; SLAMF4;

Background:

This gene encodes a cell surface receptor expressed on natural killer (NK) cells (and some T cells) that mediate non-major histocompatibility complex (MHC) restricted killing. The interaction between NK-cell and target cells via this receptor is thought to modulate NK-cell cytolytic activity. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).

Catalog #:A1993

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:51744

Isotype:IgG

Swiss Prot:Q9BZW8

Purity:Affinity purification

For research use only.