

AP2B1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:105kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human AP2B1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

i

Synonym:

ADTB2; AP105B; CLAPB1; AP2-BETA

Catalog #:A1995

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:163

Isotype:IgG

Swiss Prot:P63010

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is one of two large chain components of the assembly protein complex 2, which serves to link clathrin to receptors in coated vesicles. The encoded protein is found on the cytoplasmic face of coated vesicles in the plasma membrane. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).