

AOC3

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:84kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human AOC3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HPAO; SSAO; VAP1; VAP-1;

Background:

This gene encodes a member of the semicarbazide-sensitive amine oxidase family. Copper amine oxidases catalyze the oxidative conversion of amines to aldehydes in the presence of copper and quinone cofactor. The encoded protein is localized to the cell surface, has adhesive properties as well as monoamine oxidase activity, and may be involved in leukocyte trafficking. Alterations in levels of the encoded protein may be associated with many diseases, including diabetes mellitus. A pseudogene of this gene has been described and is located approximately 9-kb downstream on the same chromosome. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).

Catalog #:A2001

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8639

Isotype:IgG

Swiss Prot:Q16853

Purity:Affinity purification

For research use only.