

TNFSF14

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:1000 - 1:2000 IHC 1:50 - 1:200

Calculated MW:26kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TNFSF14

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LTg; TR2; CD258; HVEM; LIGHT;

Background:

The protein encoded by this gene is a member of the tumor necrosis factor (TNF) ligand family. This protein is a ligand for TNFRSF14, which is a member of the tumor necrosis factor receptor superfamily, and which is also known as a herpesvirus entry mediator (HVEM). This protein may function as a costimulatory factor for the activation of lymphoid cells and as a deterrent to infection by herpesvirus. This protein has been shown to stimulate the proliferation of T cells, and trigger apoptosis of various tumor cells. This protein is also reported to prevent tumor necrosis factor alpha mediated apoptosis in primary hepatocyte. Two alternatively spliced transcript variant encoding distinct isoforms have been reported.

To place an order, please [Click HERE](#).

Catalog #:A2002

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8740

Isotype:IgG

Swiss Prot:O43557

Purity:Affinity purification

For research use only.