

## POSTN

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**93kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human POSTN

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

PN; OSF2; OSF-2; PDLPOSTN; periostin;

*To place an order, please [Click HERE](#).*

**Catalog #:**A2003

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**10631

**Isotype:**IgG

**Swiss Prot:**Q15063

**Purity:**Affinity purification

For research use only.