

IGFBP4

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:28kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human IGFBP4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BP-4; IBP4; IGFBP-4; HT29-IGFBP;

Catalog #:A2008

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3487

Isotype:IgG

Swiss Prot:P22692

Purity:Affinity purification

For research use only.

Background:

This gene is a member of the insulin-like growth factor binding protein (IGFBP) family and encodes a protein with an IGFBP domain and a thyroglobulin type-I domain. The protein binds both insulin-like growth factors (IGFs) I and II and circulates in the plasma in both glycosylated and non-glycosylated forms. Binding of this protein prolongs the half-life of the IGFs and alters their interaction with cell surface receptors.

To place an order, please [Click HERE](#).