

EIF5A

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:16kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human EIF5A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

EIF-5A; EIF5A1; eIF5A1;

Catalog #:A2016

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1984

Isotype:IgG

Swiss Prot:P63241

Purity:Affinity purification

For research use only.

Background:

mRNA-binding protein involved in translation elongation. Has an important function at the level of Mrna turnover, probably acting downstream of decapping. Involved in actin dynamics and cell cycle progression, mRNA decay and probably in a pathway involved in stress response and maintenance of cell wall integrity. With syntenin SDCBP, functions as a regulator of p53/TP53 and p53/TP53-dependent apoptosis. Regulates also TNF-alpha-mediated apoptosis. Mediates effects of polyamines on neuronal process extension and survival. May play an important role in brain development and function, and in skeletal muscle stem cell differentiation. Also described as a cellular cofactor of human T-cell leukemia virus type I (HTLV-1) Rex protein and of human immunodeficiency virus type 1 (HIV-1) Rev protein, essential for mRNA export of retroviral transcripts.

To place an order, please [Click HERE](#).