

## MT2A

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**7kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human MT2A

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

h

**Synonym:**

MT2; MT2A;MT-II;

**Catalog #:**A2018

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**4502

**Isotype:**IgG

**Swiss Prot:**P02795

**Purity:**Affinity purification

For research use only.

**Background:**

MT2A (metallothionein 2A) is a protein-coding gene. Diseases associated with MT2A include cerebral lymphoma. GO annotations related to this gene include drug binding. An important paralog of this gene is MT1G.

*To place an order, please [Click HERE](#).*