

MT2A

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 7kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MT2A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

h

Synonym:

MT2; MT2A; MT-II;

Catalog #: A2018

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4502

Isotype: IgG

Swiss Prot: P02795

Purity: Affinity purification

For research use only.

Background:

MT2A (metallothionein 2A) is a protein-coding gene. Diseases associated with MT2A include cerebral lymphoma. GO annotations related to this gene include drug binding. An important paralog of this gene is MT1G.

To place an order, please [Click HERE](#).