

CXCL5

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:12kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CXCL5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SCYB5; ENA-78

Catalog #:A2020

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6374

Isotype:IgG

Swiss Prot:P42830

Purity:Affinity purification

For research use only.

Background:

This gene encodes a protein that is a member of the CXC subfamily of chemokines. Chemokines, which recruit and activate leukocytes, are classified by function (inflammatory or homeostatic) or by structure. This protein is proposed to bind the G-protein coupled receptor chemokine (C-X-C motif) receptor 2 to recruit neutrophils, to promote angiogenesis and to remodel connective tissues. This protein is thought to play a role in cancer cell proliferation, migration, and invasion.

To place an order, please [Click HERE](#).