

PRDX6

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:25kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PRDX6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

1-Cys; AOP2; KIAA0106; MGC46173; NSGPx; PRX; aiPLA2; p29;

Background:

Prdx6(Peroxiredoxin-6) , a unique member of the peroxiredoxin family of antioxidants, is highly expressed in liver and protects cells from oxidative damage by reducing H₂O₂ and various lipid Peroxides(PMID: 17382207).It is also named as AOP2, KIAA0106. It can form a dimer(PMID:20500660).PRDX6 is expressed in all major organs, with a particularly high level in lung(PMID:15890616).This antibody is specific to PRDX6.

To place an order, please [Click HERE](#).

Catalog #:A2031

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9588

Isotype:IgG

Swiss Prot:P30041

Purity:Affinity purification

For research use only.