

CD83

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:23kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CD83

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BL11; HB15;

Background:

The protein encoded by this gene is a single-pass type I membrane protein and member of the immunoglobulin superfamily of receptors. The encoded protein may be involved in the regulation of antigen presentation. A soluble form of this protein can bind to dendritic cells and inhibit their maturation. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]

To place an order, please [Click HERE](#).

Catalog #:A2040

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9308

Isotype:IgG

Swiss Prot:Q01151

Purity:Affinity purification

For research use only.