

CD69

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:20kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CD69

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AIM; EA1; MLR-3; CLEC2C; GP32/28; BL-AC/P26;

Background:

This gene encodes a member of the calcium dependent lectin superfamily of type II transmembrane receptors. Expression of the encoded protein is induced upon activation of T lymphocytes, and may play a role in proliferation. Furthermore, the protein may act to transmit signals in natural killer cells and platelets. [provided by RefSeq, Aug 2011]

To place an order, please [Click HERE](#).

Catalog #:A2045

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:969

Isotype:IgG

Swiss Prot:Q07108

Purity:Affinity purification

For research use only.