

H2AFX

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF IP

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200 IF 1:20 - 1:100 IP 1:20 - 1:100

Calculated MW:15kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human H2AFX

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

H2A.X; H2A/X; H2AX;

Background:

H2AFX (H2A histone family, member X) is a protein-coding gene. Diseases associated with H2AFX include vulvar intraepithelial neoplasia. GO annotations related to this gene include histone binding and protein heterodimerization activity. An important paralog of this gene is HIST3H2A.

To place an order, please [Click HERE](#).

Catalog #:A2082

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3014

Isotype:IgG

Swiss Prot:P16104

Purity:Affinity purification

For research use only.