

REG3A

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:200

Calculated MW:19kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human REG3A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HIP; INGAP; PAP; PAP-H; PAP1; PBCGF; REG-III; REG3;

Background:

This gene encodes a pancreatic secretory protein that may be involved in cell proliferation or differentiation. It has similarity to the C-type lectin superfamily. The enhanced expression of this gene is observed during pancreatic inflammation and liver carcinogenesis. Alternate splicing results in multiple transcript variants that encode the same protein.[provided by RefSeq, Dec 2010]

To place an order, please [Click HERE](#).

Catalog #:A2101

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5068

Isotype:IgG

Swiss Prot:Q06141

Purity:Affinity purification

For research use only.