

NUMBL

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:65kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human NUMBL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

t

Synonym:

CAG3A; CTG3a; NBL; NUMB-R; NUMBLIKE; NUMBR; TNRC23;

Catalog #:A2140

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9253

Isotype:IgG

Swiss Prot:Q9Y6R0

Purity:Affinity purification

For research use only.

Background:

Numb-like (Numbl or Nbl) binds Notch and is an intracellular modulator of Notch signaling pathway. Numbl contains a phosphotyrosine binding (PTB) domain that could bind the Y682 motif present in the cytoplasmic tail of APP. Deduced molecular weight of Numbl is 72 kDa.

To place an order, please [Click HERE](#).