

## NUMBL

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**65kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human NUMBL

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

t

**Synonym:**

CAG3A; CTG3a; NBL; NUMB-R; NUMBLIKE; NUMBR; TNRC23;

**Catalog #:**A2140

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**9253

**Isotype:**IgG

**Swiss Prot:**Q9Y6R0

**Purity:**Affinity purification

For research use only.

**Background:**

Numb-like (Numbl or Nbl) binds Notch and is an intracellular modulator of Notch signaling pathway. Numbl contains a phosphotyrosine binding (PTB) domain that could bind the Y682 motif present in the cytoplasmic tail of APP. Deduced molecular weight of Numbl is 72 kDa.

*To place an order, please [Click HERE](#).*